

INTRODUCER KITS



Stay on the right track.

Introducer Kits by Typenex® Medical are designed to assist with smooth vascular access. Featuring a needle, guidewire, sheath, and dilator, this all-in-one system offers ease and efficiency. With a variety of sizes and styles available, you can choose the solution that best suits each procedure.



FEATURES AND BENEFITS

Efficient Design — All-in-one system consolidates inventory and streamlines workflows.

Multiple Sizes — Various options to accommodate different procedures, preferences, and patients.

Tailored Components — Echogenic needles and marked specialty sheaths available for added ease and convenience.

Hemostasis Valve — Optional HV valve helps minimize blood loss during vasculature introduction.

STANDARD INTRODUCER KITS

Part Number	Description	Needle	Guidewire	HV Valve	Units/Case
MS0203	4 Fr x 12 cm Elite HV Micro-Introducer Kit	Standard	Stainless Steel	Yes	10
MS0303	4 Fr x 24 cm Elite HV Micro-Introducer Kit	Standard	Stainless Steel	Yes	5
MS0403	4 Fr x 45 cm Elite HV Micro-Introducer Kit	Standard	Stainless Steel	Yes	10
MS0511	4 Fr x 70 cm Elite HV Micro-Introducer Kit	Standard	Nitinol	Yes	10
MS0603	5 Fr x 10 cm Coax Micro-Introducer Kit	Standard	Stainless Steel	No	10

ECHOGENIC INTRODUCER KITS

Part Number	Description	Needle	Guidewire	HV Valve	Units/Case
MS0236	4 Fr x 12 cm Elite HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	10
MS0436	4 Fr x 45 cm Elite HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	10
MS0636	5 Fr x 10 cm Coax Micro-Introducer Kit	Echogenic	Nitinol	No	10
VASNIS02A	6 Fr x 7 cm HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	5
VASNIS03A	7 Fr x 7 cm HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	5
VASMI4F01	4 Fr x 10 cm Micro-Puncture Kit	Echogenic	Nitinol	No	5

STIFFEN INTRODUCER KITS

Part Number	Description	Needle	Guidewire	HV Valve	Units/Case
MS1235S	4 Fr x 30 cm Coax Stiffen Micro-Introducer Kit	Standard	Nitinol	No	10
MS1335S	4 Fr x 40 cm Coax Stiffen Micro-Introducer Kit	Standard	Nitinol	No	5
MS1835S	4 Fr x 60 cm Coax Stiffen Micro-Introducer Kit	Standard	Nitinol	No	10

RADIOFREQUENCY KITS

Part Number	Description	Needle	Guidewire	HV Valve	Units/Case
MS1512	6 Fr x 12 cm HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	10
MS2012	6 Fr x 7 cm Coax Micro-Introducer Kit	Echogenic	Nitinol	No	10
MS2212	7 Fr x 7 cm Coax Micro-Introducer Kit	Echogenic	Nitinol	No	10

Manufactured for Typenex® Medical.

Visit www.typenex.com or call 866-897-3639 to learn more.