

Fluid Fighter[®] Suction Canisters and Liners

Clean. Compact. Compatible.



INTRODUCTION

Typenex® Medical is proud to offer a comprehensive portfolio of reusable suction canisters, disposable flexible liners, and auxiliary accessories for streamlined fluid waste retention and disposal in clinical settings. Designed to support standardization and minimize medical waste, this sturdy yet lightweight system is available in a variety of configurations to suit your specific fluid management needs.

FEATURES AND BENEFITS



Fluid Fighter Suction Canisters

Shatter-Resistant — Reusable canisters feature an integrated one-piece design; backed by a 20-year guarantee.

Compatible — Adapters and custom solutions offer seamless integration and system compatibility with Cardinal Health, Medline, Bemis, and DeRoyal.¹

Increased Visibility — Approximate measurement graduations are clearly marked in black text.



Fluid Fighter Suction Liners

Compact — Optimal packaging and minimal assembly required, improving turnover times and conserving valuable storage space.

Self-Seating — Suction liners seat themselves in about 8 seconds, 50% faster than competing brands.

Closed System — Built-in handles allow for quick and easy removal after use, preventing contact with fluid waste.

Durable — Liners have been drop tested² and stretch up to 270% of initial dimension.

Infection Prevention — Integrated hydrophobic filter to prevent backflow of fluids and bacterial/viral transfer to vacuum source.³

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE	QUANTITY
SC0101	Fluid Fighter Suction Liner	1000 cc	100/Case
SC0102	Fluid Fighter Suction Liner	1500 cc	100/Case
SC0103	Fluid Fighter Suction Liner	2000 cc	100/Case
SC0104	Fluid Fighter Suction Liner	3000 cc	70/Case
SC0201	Fluid Fighter Suction Canister	1000 cc	15/Case
SC0202	Fluid Fighter Suction Canister	1500 cc	15/Case
SC0203	Fluid Fighter Suction Canister	2000 cc	12/Case
SC0204	Fluid Fighter Suction Canister	3000 cc	9/Case
SC0405CVL	Vacuum Control L-Valve	—	50/Bag
SC0405CVT	Vacuum Control T-Valve	—	50/Bag
SC0405XVTC	Vacuum T-Connector	—	50/Bag
SC0405ET	Extension Tube, Female Connector	0.23 in x 5.9 ft	1/Each
SC0405LT	Suction Canister Trolley	—	1/Case
SC0405RH	Rail Hanger	—	1/Each
SC0405SCT	Serial Connection Tube	—	100/Case
SC0405ST	Silicone Tube	100 m (7.2 mm I.D.)	1/Case
SC0405WH	Wall Hanger Set	—	1/Each
SC0405WHA	Wall Hanger Adapter (Compatible with Cardinal, Medline, Bemis)	—	1/ Each
PER05TD3669	Wall Hanger Adapter (Compatible with DeRoyal)	—	1/ Each
PER05TD8770	Serial Connection Connector	—	80/Case
PER05TD0046	Suction Liner Canister Vacuum Connector	—	50/Case
PER05TD1618	Solidifying Agent, Dissolvable Pouch	—	150/Bag

Solidifiers are compatible with 1000 cc, 1500 cc, 2000 cc, and 3000 cc Fluid Fighter® Suction Liners.

¹Compatibility verified on Cardinal Health roll stands and part number 65652-145.
Compatibility verified on Medline part numbers DYNDSCWL, DYNDSCC, and ORBASE.
Compatibility verified on Bemis Health Care part number 530510.
Compatibility verified on DeRoyal Industries part numbers 71-5017, 71-1116, 71-5020, and 71-5045

²Typenex® Medical's Fluid Fighter® Suction Liners are drop tested according to ASTM D5276-98.
³Bacterial filtration efficiency of >99.85% and viral filtration efficiency of >97.5% validated by Pergo Medikal.

All products and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only.

FOR MORE INFORMATION

Visit www.typenex.com or call
866-897-3639 to learn more.



300 East Randolph Street
Suite 40.100
Chicago, IL 60601