



Sample a blood segment with ease.

Typesafe® is a segment piercing device that allows laboratory technicians to sample a blood segment without the need for scissors. This single-use solution is designed to help support safety and efficiency during pre-transfusion testing.

Typesafe® | Segment Piercing Device



Features and Benefits

Efficient — Dispenses blood from a blood segment into a collection tube without the use of scissors.

Convenient — Easy to handle with gloved hands; not made with natural rubber latex.

Stable — Fits securely onto 10 mm and 12 mm tubes; accommodates most standard blood segment sizes.



Ordering Information

Part Number	Description	Units/Case
TS500	Typesafe Segment Piercing Device	500
4R4750	Typesafe Segment Piercing Device	1000

For more information, visit us online at www.typenex.com or call us at 866-897-3639.