



Smooth streaking, smooth savings.

Typenex® Medical offers a comprehensive selection of disposable inoculation loops and needles to support a wide range of specimen streaking applications and preferences. Choose from various configurations and packaging styles to best suit your process while securing substantial cost savings.



Inoculation Loops and Needles

Inoculation Loops and Needles



Features and Benefits

Versatile Options — Choose soft, flexible loops for gentle malleability or rigid formats for firm control.

Thoughtful Design — Hexagonal shaft provides improved grip and orientation; smooth loop surfaces are free of burrs and flashing.

Convenient Packaging — Available in different packaging styles such as easy-to-tear perforated bags and resealable bags.

Ready to Use — Sterile, calibrated,* color-coded, and offered in various configurations for seamless integration with existing workflows.



Perforated Bag

Ordering Information

Part Number	Description	Format	Size	Color	Packaging Type	Units/Case
LP0102	Inoculation Loops	Rigid	1 μ L	Orange	Resealable	1000
LP0202	Inoculation Loops	Rigid	10 μ L	Dark Blue	Resealable	1000
LP0103	Inoculation Loops	Flexible	1 μ L	Yellow	Resealable	1000
LP0203	Inoculation Loops	Flexible	10 μ L	Green	Resealable	1000
LP0302	Inoculation Needles	—	—	Red	Resealable	1000
LP0102PF	Inoculation Loops	Rigid	1 μ L	Orange	Perforated	1000
LP0202PF	Inoculation Loops	Rigid	10 μ L	Dark Blue	Perforated	1000
LP0103PF	Inoculation Loops	Flexible	1 μ L	Yellow	Perforated	1000
LP0203PF	Inoculation Loops	Flexible	10 μ L	Green	Perforated	1000
LP0302PF	Inoculation Needles	—	—	Red	Perforated	1000

*Calibrated using the Evans Blue Dye Method.

Visit www.typenex.com or call 866-897-3639 to learn more.