



# Affordable and universally compatible. How rad.

Typenex® Medical is proud to introduce Rad-Control Universal, an innovative indicator tag that serves to verify the irradiation of blood products by both gamma and x-ray sources. This affordable, universally compatible solution is designed to help consolidate resources and seamlessly accommodate existing workflows.

# Rad-Control Universal | Irradiation Tag



## Features and Benefits

**Consolidation** — Compatible with both gamma and x-ray sources.

**Visual Verification** — Serves to verify the irradiation of blood products by providing visual confirmation of exposure.

**Customization** — Add facility logo, custom prompts, and optional logbook stickers to best suit existing workflows.

**Cost-Savings** — Competitively priced to help secure substantial savings.

## Ordering Information

Part Number	Description	Irradiation Type	Units/Case
RT0301	Rad-Control   Universal	25 Gy, Universal	200
RT0302	Rad-Control   Universal	15 Gy, Universal	200

Rad-Control Universal should not be considered as a radiation dosimeter

Visit [www.typenex.com](http://www.typenex.com) or call 866-897-3639 to learn more.