

INTRODUCER KITS



Stay on the right track.

Introducer Kits by Typenex® Medical are designed to assist with smooth vascular access. Featuring a needle, guidewire, sheath, and dilator, this all-in-one system offers ease and efficiency. With a variety of sizes and styles available, you can choose the solution that best suits each procedure.



FEATURES AND BENEFITS

Efficient Design — All-in-one system consolidates inventory and streamlines workflows.

Multiple Sizes — Various options to accommodate different procedures, preferences, and patients.

Tailored Components — Echogenic needles and marked specialty sheaths available for added ease and convenience.

Hemostasis Valve — Optional HV valve helps minimize blood loss during vasculature introduction.

4 FR INTRODUCER KITS

Part Number	Description	Size	Needle	Guidewire	HV Valve	Units/Case
MS0203	Elite HV Micro-Introducer Kit	4 Fr x 12 cm	Standard	Stainless Steel	Yes	10
MS0236	Elite HV Micro-Introducer Kit	4 Fr x 12 cm	Echogenic	Nitinol	Yes	10
MS0403	Elite HV Micro-Introducer Kit	4 Fr x 45 cm	Standard	Stainless Steel	Yes	10
MS1235S	Coax Stiffen Micro-Introducer Kit	4 Fr x 30 cm	Standard	Nitinol	No	10
MS1335S	Coax Stiffen Micro-Introducer Kit	4 Fr x 40 cm	Standard	Nitinol	No	5
MS1835S	Coax Stiffen Micro-Introducer Kit	4 Fr x 60 cm	Standard	Nitinol	No	10
VASRS041118S-II	HV Introducer Kit	4 Fr x 11 cm	Standard	Stainless Steel	Yes	5

5 FR INTRODUCER KITS

Part Number	Description	Size	Needle	Guidewire	HV Valve	Units/Case
MS0603	Coax Micro-Introducer Kit	5 Fr x 10 cm	Standard	Stainless Steel	No	10
MS0636	Coax Micro-Introducer Kit	5 Fr x 10 cm	Echogenic	Nitinol	No	10
VASRS051118S-II	HV Introducer Kit	5 Fr x 11 cm	Standard	Stainless Steel	Yes	5

6 FR INTRODUCER KITS

Part Number	Description	Size	Needle	Guidewire	HV Valve	Units/Case
MS1512	HV Micro-Introducer Kit	6 Fr x 12 cm	Echogenic	Stainless Steel	Yes	10
VASRS060718S-II	HV Introducer Kit	6 Fr x 7 cm	Standard	Stainless Steel	Yes	5


7 FR INTRODUCER KITS


Part Number	Description	Size	Needle	Guidewire	HV Valve	Units/Case
MS2212	Coax Micro-Introducer Kit	7 Fr x 7 cm	Echogenic	Nitinol	No	10
VASRS070718S-II	HV Introducer Kit	7 Fr x 7 cm	Standard	Stainless Steel	Yes	5
VASRS071118S-II	HV Introducer Kit	7 Fr x 11 cm	Standard	Stainless Steel	Yes	5

MS-series kits are manufactured for Typenex® Medical.

VASR-series kits are manufactured by Vascular Access Direct and distributed by Typenex® Medical.

Visit www.typenex.com or call 866-897-3639 to learn more.

 300 East Randolph Street
Suite 40.100
Chicago, IL 60601

 866-897-3639
866-897-3630 (fax)

 sales@typenex.com

 www.typenex.com