



The clear choice in bulb syringes.

Visibulb® allows the user to see the contents suctioned from ear, oral, and nasal passageways. Unlike traditional bulb syringes made of opaque materials, Visibulb's clear design offers transparency when suctioning the patient's orifice.

Visibulb® | 2 oz



Features and Benefits

- Clear 2 oz bulb syringe designed to remove fluids from ear, oral, and nasal passageways with a simple squeeze
- Flat bottom allows the bulb syringe to stand without tipping
- Not made with natural rubber latex
- Sterile, single-use

Instructions

1. Squeeze the bulb of the syringe, then place it gently inside the patient's orifice.
2. When the pressure has been completely released from the bulb, remove it from the patient's orifice.

Note: You should see fluids being pulled into the bulb of the syringe.

3. Repeat the prior steps until the patient's orifice is clear of fluids.

Ordering Information

Part Number	Description	Units/Case
CEB02	Visibulb 2 oz Ear and Ulcer Syringe	50

Visibulb® is not intended to be used as a diagnostic tool and should not be used to determine course of action for patient care.

Visit www.typenex.com or call 866-897-3639 to learn more.