



Smooth streaking, smooth savings.


Typenex® Medical offers a comprehensive selection of disposable inoculation loops and needles to support a wide range of specimen streaking applications and preferences. Choose from various configurations and packaging styles to best suit your process while securing substantial cost savings.




Perforated



Resealable



Tube



Features and Benefits

Versatile Options — Choose soft, flexible loops for gentle malleability or rigid formats for firm control.

Thoughtful Design — Hexagonal shaft provides improved grip and orientation; smooth loop surfaces are free of burrs and flashing.

Convenient Packaging — Available in different packaging styles such as resealable options, easy-to-tear perforated bags, and free-standing plastic tubes.

Ready to Use — Sterile, calibrated,* color-coded, and offered in various configurations for seamless integration with existing workflows.

Ordering Information

Part Number	Description	Packaging Type	Format	Size	Color	Units/Case
LP0103PF	Inoculation Loops	Perforated	Flexible	1 μ L	Yellow	1000
LP0203PF	Inoculation Loops	Perforated	Flexible	10 μ L	Green	1000
LP0102PF	Inoculation Loops	Perforated	Rigid	1 μ L	Orange	1000
LP0202PF	Inoculation Loops	Perforated	Rigid	10 μ L	Dark Blue	1000
LP0302PF	Inoculation Needles	Perforated	—	—	Red	1000
LP0103	Inoculation Loops	Resealable	Flexible	1 μ L	Yellow	1000
LP0203	Inoculation Loops	Resealable	Flexible	10 μ L	Green	1000
LP0102	Inoculation Loops	Resealable	Rigid	1 μ L	Orange	1000
LP0202	Inoculation Loops	Resealable	Rigid	10 μ L	Dark Blue	1000
LP0302	Inoculation Needles	Resealable	—	—	Red	1000
LP0103TB	Inoculation Loops	Tube	Flexible	1 μ L	Yellow	1000
LP0203TB	Inoculation Loops	Tube	Flexible	10 μ L	Green	1000
LP0102TB	Inoculation Loops	Tube	Rigid	1 μ L	Orange	1000
LP0202TB	Inoculation Loops	Tube	Rigid	10 μ L	Dark Blue	1000
LP0302TB	Inoculation Needles	Tube	—	—	Red	1000

*Calibrated using the Evans Blue Dye Method.

Visit www.typenex.com or call 866-897-3639 to learn more.

300 East Randolph Street
Suite 40.100
Chicago, IL 60601

866-897-3639
866-897-3630 (fax)

sales@typenex.com

www.typenex.com