




Take a stand against floor fluid buildup.

The Fluid Fighter® Suction Mat by Typenex® Medical is thoughtfully engineered to capture fluid from the surface while suction is drawn through a system of channels within the mat. This low-profile floor fluid management solution is designed with high-quality materials and can be securely adhered to your operating room floor.

Fluid Fighter® | Suction Mat

TYPENEX
MEDICAL

Fluid Fighter® | Suction Mat



Features and Benefits

Dual Action — Absorbent top layer works simultaneously with an internal suction system to better manage fluid buildup.

Low-Profile — Four adhesive patches grip to the floor, keeping the mat flat and bulk-free.


Quality Design — Size and features are optimized to provide the best coverage and suction.


Efficient — Easy cleanup and disposal allows for a quick turnover between procedures.

Ordering Information

Part Number	Description	Size	Units/Case
FM0202QD	Fluid Fighter I Suction Mat	47" x 35"	12

Visit www.typhenex.com or call 866-897-3639 to learn more.

 300 East Randolph Street
Suite 40.100
Chicago, IL 60601

 866-897-3639
866-897-3630 (fax)

 sales@typhenex.com

 www.typhenex.com