




Your operating room floor's first line of defense.

The Fluid Fighter® Suction Strip is a sleek, low-profile device that quietly and continuously suctions fluid from the operating room floor. Designed with strategically placed suction ports and absorbent padding underneath, the Fluid Fighter Suction Strip can be easily applied and repositioned to suction, absorb, and contain floor fluid buildup — helping to improve room turnover time.

Fluid Fighter[®] | Suction Strip



Features and Benefits

Quiet — Absorbent material on the bottom of the device helps mitigate noise.

Powerful — Continuously suctions floor fluid waste throughout fluid-intensive procedures.


Compatible — Pre-attached 12-ft tubing suitable for all standard operating room suction systems.




Ordering Information

Part Number	Description	Size	Units/Case
FM0701	Fluid Fighter Suction Strip	31 in x 2 in	10

Visit www.tylenex.com or call 866-897-3639 to learn more.

 300 East Randolph Street
Suite 40.100
Chicago, IL 60601

 866-897-3639
866-897-3630 (fax)

 sales@tylenex.com

 www.tylenex.com