



# Shhhh...floor suction can be distracting.

Introducing the Big Foot™ Quiet by Typenex® Medical, our quieter solution to floor suction. This enhanced version of our Big Foot™ Suction Disc has an absorbent layer on the underside of the disc. This feature mitigates unwanted noise and helps effectively vacuum fluid off the floor. A dry, clean operating room is an effective one.

Big Foot™ Quiet | Suction Disc

TYPENEX   
MEDICAL

# Big Foot™ Quiet | Suction Disc



## Features and Benefits

**Quiet** — Bottom padding helps eliminate disruptive noise that can interfere with vital communication between clinicians.

**Suction and Absorption** — Absorbent padding on the underside of the disc acts as an acquisition layer, helping to catch and contain fluid before suctioning.

**16 Suction Ports** — Ports are strategically distributed across the device to maximize fluid removal.

**Sleek Design** — Side port connection to tubing allows for ease of maneuverability across the floor to address areas of high fluid volume.


**Compatible** — Suitable for all standard OR suction systems through a simple attachment of tubing.

## Ordering Information

Part Number	Description	Color	Units/Case
FM0301Q	Big Foot Quiet I Suction Disc	Black	40

Visit [www.typenex.com](http://www.typenex.com) or call 866-897-3639 to learn more.

 300 East Randolph Street  
Suite 40.100  
Chicago, IL 60601

 866-897-3639  
866-897-3630 (fax)

 [sales@typenex.com](mailto:sales@typenex.com)

 [www.typenex.com](http://www.typenex.com)