



Locked in for  
patient verification.

Designed to help reduce blood transfusion errors associated with patient identification, FinalCheck® is an easy-to-use mechanical barrier system that can only be opened using a code that is placed on the patient's admissions band.

FinalCheck® | Recipient Verification System

TYPENEX  
MEDICAL



## Features and Benefits

**Patient Safety** — With the FinalCheck Recipient Verification System, you can ensure that pre-transfusion patient ID checks are done at the bedside by putting this extra layer of safety in place.

**Efficient** — This physical barrier to the blood component can only be opened by entering the code that is applied to the patient's admissions band, ensuring the right blood goes to the right patient.

**Ease of Use** — The FinalCheck lock is easy to set and open using the four-digit code on the patient label set.

## Ordering Information

Part Number	Color	Detail	Units/Case
<b>LOCK</b>			
C3000	—	—	100
<b>LABELS</b>			
C3100	Red	1 Specimen Tube Sticker	3000
C3103	Blue	1 Specimen Tube Sticker	3000
C3113	Blue	2 Specimen Tube Stickers	3000
<b>BAGS</b>			
C3201	Clear	Biohazard Symbol	100
<b>BANDS</b>			
C4678	Red	Additional Stickers	3000

Visit [www.typenex.com](http://www.typenex.com) or call 866-897-3639 to learn more.