



# Help warm up the heel stick process.

The Infant Heel Warmer from Typenex® Medical is designed to deliver the right amount of warmth to increase blood circulation in an infant's heel during blood sampling. This single-use solution provides even heat distribution and consistent warming to help caregivers successfully collect a heel stick blood draw.

## Infant Heel Warmer

# Infant Heel Warmer



## Features and Benefits

**Effective** — Designed to conform to the infant's foot shape, keeping warmth where needed.

**Convenient** — Features a repositionable adhesive strip to keep the heel warmer in place.

**Easy to Use** — Simple disk activation offers convenient, streamlined preparation.

**Reliable** — Delivers 103 °F of heat therapy when activated at room temperature; will not exceed 105 °F.


**Safe** — Made with non-toxic, patented Insta-Gel technology.

## Ordering Information

Part Number	Description	Color	Units/Case
HWA94245	Infant Heel Warmer	Clear	100

Visit [www.typenex.com](http://www.typenex.com) or call 866-897-3639 to learn more.

 300 East Randolph Street  
Suite 40.100  
Chicago, IL 60601

 866-897-3639  
866-897-3630 (fax)

 [sales@typenex.com](mailto:sales@typenex.com)

 [www.typenex.com](http://www.typenex.com)