

MICROTOME BLADES

DuraEdge®

- Manufactured using a proprietary process for hardened stainless steel, ensuring a finely honed blade with a flawless edge
- Capable of producing ultra-thin sections of tissue
- Uniform PTFE coating
- 35° blade angle



Part Number	Description	Profile	Dimensions	Dispenser Color	Units/Case
SMCS7223	DuraEdge Microtome Blades	Low	77 x 8 mm	Green	500
SMCS7310	DuraEdge Microtome Blades	High	77 x 14 mm	Red	500

Encore™

- Designed with enhanced sharpness, easily cutting up to three blocks without overlapping
- Produces consistent sectioning performance
- Uniform PTFE coating
- 35° blade angle



Part Number	Description	Profile	Dimensions	Dispenser Color	Units/Case
SMCS7380	Encore Microtome Blades	Low	80 x 8 mm	Pink, Green	500
SMCS7390	Encore Microtome Blades	High	80 x 14 mm	Pink, Green	500

Sapphire®

- Engineered with optimum thickness and precise edge geometry
- Patented Sapphire coating produces quality ribbons in fewer cuts, extending the blade life and saving time
- No-hole design provides maximum strength and precision
- 35° blade angle



Part Number	Description	Profile	Dimensions	Dispenser Color	Units/Case
SMCS8000	Sapphire Microtome Blades	Low	80 x 8 mm	Blue	500
SMCS8100	Sapphire Microtome Blades	High	77 x 14 mm	Blue	500

Visit www.typenex.com or call 866-897-3639 to learn more.