


## An innovative approach to irradiation verification.

Typenex<sup>®</sup> Medical's Rad-Control Irradiation Tags serve to verify the irradiation of blood products while supporting and streamlining workflows. With three distinct options to choose from, you can discover the solution that best suits your facility's needs.



## Features and Benefits

Rad-Control Irradiation Tags are all equipped with verification technology, easy-peel adhesives, customizable fields, and optional logbook stickers. With no need for refrigeration or manual data entry, our indicator tags can help conserve storage space and reduce the risk of transcription errors.



Rad-Control Standard after irradiation

### Rad-Control Standard

- Transitions from red to black, indicating whether the blood product has been irradiated.
- Available in configurations for either gamma or x-ray irradiation sources.



Rad-Control Barcode after Irradiation

### Rad-Control Universal

- Transitions from red to black, indicating whether the blood product has been irradiated.
- Compatible with both gamma and x-ray irradiation sources.



Rad-Control Scale after irradiation

### Rad-Control Scale

- Provides visual confirmation that the absorbed dose has not exceeded the maximum permissible level.
- Semi-quantitative, color-changing design approximates the amount of radiation dose delivered.

## Ordering Information

| Number | Description                                          | Units/Case |
|--------|------------------------------------------------------|------------|
| RT0101 | Rad-Control Irradiation Tag   Standard, 25 Gy, Gamma | 200        |
| RT0102 | Rad-Control Irradiation Tag   Standard, 15 Gy, Gamma | 200        |
| RT0201 | Rad-Control Irradiation Tag   Standard, 25 Gy, X-ray | 200        |
| RT0202 | Rad-Control Irradiation Tag   Standard, 15 Gy, X-ray | 200        |
| RT0301 | Rad-Control Irradiation Tag   Universal, 25 Gy       | 200        |
| RT0302 | Rad-Control Irradiation Tag   Universal, 15 Gy       | 200        |
| RT0303 | Rad-Control Irradiation Tag   Scale, 15–50 Gy        | 200        |
| RT0304 | Rad-Control Irradiation Tag   Scale, 25–50 Gy        | 200        |

Rad-Control Irradiation Tags should not be considered as radiation dosimeters.

Visit [www.typenex.com](http://www.typenex.com) or call 866-897-3639 to learn more.