

# Beyond Wireless, True Mobility



#### True **Mobility**

Freedom to scan anytime, anywhere.



#### **Minimal Learning Curve**

Buy today, scan tomorrow. Intuitive UI eliminates the need for training, just open and scan.



#### **Instant Setup and No Downtime**

Start scanning the moment the app opens. Always up to date-no OS or driver headaches.



#### **Engage Patients. Empower Teams.**

Show scans, explain clearly, boost treatment acceptance. Share cases and tasks in real time-work as one, anywhere.



#### **Medit Link Ecosystem**

From scan to consultation to order-stay connected every step with Medit Link & Medit Link Express.





#### **Clinicians**



# Focus on Clinical Freedom and Comfort

Scan anywhere for optimal ergonomics.

Position the scanner for comfort, reduce physical strain.

#### **Patients**



**Clarity & Confidence** 

Show and explain with real scan visuals. Patients understand instantly.

#### **Dental Team**




**Minimal Learning Curve** 

No specialists required. Any staff can easily start scanning right away.

# **True Mobility, Your Way**

Enabling Your Vision and Next-Level Experiences with i900 Mobility.



• Use one device across multiple chairs in busy practices



• Community center pop-up clinics and field screening events



• Deliver live scan demos at conferences, seminars, or off-site events



**Remote Work & After-Hours Access** 

• Approve or place orders on the go — even from the car

# MEDIT i 900 Mobility







# With Medit Link Express on iPad, Build Your Digital Dentistry Journey

#### 1st iPad-Native Standalone App

A dedicated app built exclusively for iPad, optimized for mobile intraoral scanning workflows.



#### **Optimized for Scan & Send**

Focused workflow from patient registration to scanning and sending, enabling fast clinical turnaround.



#### **Simple and Intuitive User Experience**

Streamlined interface with only essential functions, enabling quick adoption with minimal learning curve.



#### **Fast Setup, Minimal Maintenance**

Short startup time and automatic updates reduce IT overhead and allow seamless usage.



#### **Seamless Medit Link Synchronization**

All case data is automatically synced with Medit Link, enabling seamless handoff between iPad and desktop.



#### **Essential Features Built-In**

Easily complete core actions like account sign-up and partnership requests directly within Medit Link Express.



# A Native iPad Experience:

SCAN.
ENGAGE.
SEND.



- \* The iPad is not bundled with the Medit i900 Mobility; customers are advised to purchase it separately if needed.
- \* Recommended iPad model: M4 (11 inch or 13 inch)

# Medit i900 Mobility

#### Powering Accuracy with Innovation

#### 3rd-Gen Optics, Sharper Than Ever

Advanced 10-bit imaging for ultra-realistic scan results.

- Enhanced clarity with higher bit depth
- True-to-life color reproduction
- High-resolution texture recognition



Medit 10-bit imaging technology

Comparative Product

#### Wider View, Faster Workflows

18×15mm FOV and 30mm depth for broader, deeper capture.


- Capture more per pass
- · Ideal for complex or posterior regions
- · Reduces scan time

# Medit i900 Mobility

Comparative Product

#### **Sharper, Smarter Metal Scanning**

Built-in Smart filters for accurate scans on reflective surfaces.



- Metal-optimized filtering
- Reduced noise, clearer surfaces
- Reliable even in complex prosthetic

#### Glossy Mode

Revolutionary rendering for Ultra-Realistic.



Realistic 3D visualization

**All Medit scanners are designed to evolve -** regular software updates unlock new features and performance enhancements, ensuring your device keeps getting better over time.

### Smarter Hardware. Greater Possibilities.

Wifi 6E

#### Two Ways to Connect, **One Seamless Experience**

Hub or Wi-Fi? i900M fits your workflow.

- Choose traditional hub or Direct Wi-Fi Connection
- iPad supports Direct Connection only
- Flexible setup across devices and clinics



#### **Smart Battery Management, All-Day Performance**

Efficiency, protection, and power—built in.

- Up to 1 hour of continuous scanning (\*It may vary depending on the user's environment.)
- · Smart charging reduces heat and extends battery life
- Real-time battery alerts in Medit Link apps
- Auto-protection from overheat or voltage issues




Battery Charger LED

Real-Time Alerts Example

#### Fast, Wireless, Worry-Free

Supports the latest Wi-Fi 6E standard.



	Wi-Fi 5	Wi-Fi 6	Wi-Fi 6E
Frequency Bands	5GHz	2.4GHz +5GHz	2.4GHz +5GHz + NEW 6GHz
Maximum Speed	Up to 3.5Gbps	Up to 9.6Gbps	Up to 9.6Gbps with more bandwidth & less interference


#### Intuitive Scan & Control Buttons

Effortless operation at your fingertips.

- Start/stop with tactile silicone buttons
- Customizable hotkeys via Medit Link
- · Pan, zoom, rotate—without touching the screen

**Scan Button** 


**Control Button** 



# **Medit Link Express**


Effortless Mobile Scanning

A standalone app for iPad. Scan and send instantly, anytime, anywhere.




- Start scanning with Express Scan Mode or with Treatment Scan Mode depending on your workflow preference.
- Simple workflows tailored to your needs.

Scan



- Quick and intuitive scanning with a streamlined UI.
- Essential tools only—designed for speed and simplicity.
- Real-time screen sharing enhances patient consultations.

Send



- After scanning, easily send orders with just a few taps.
- Add attachments and notes to include key case details.

#### A Paradigm Shift in Dental Workflow

Every Patient Deserves Instant Clarity.

# **Express Scan Mode**

Fast Scans, Instant Conversations

Start scanning instantly.

No patient registration required.

This mode is faster, helping patients see and understand their oral condition in real time.


For visual reference use only. Not intended for prosthetic fabrication





Ready Today, Smarter Tomorrow.

Ready out of the box with automatic updates and ongoing improvements.





Rebuilt for Blazing Speed & Performance



Show 3D scans instantly on iPad and build transparent trust.




Increase treatment acceptance through interactive communication

Show More, Earn Trust, Deliver Faster.

2

# **True Mobility, Your Way**

Choose the Version That Fits Your Needs—Light or Full, We've Got You Covered.



#### **Medit Link Express - Seamless Medit Link Synchronization**

#### **One Account. Fully Connected.**

Use the same Medit account to access both Medit Link and Express. Your data stays automatically synced across platforms.

#### **Built for Workflow Flexibility**

Start with quick scans and chairside consultations in Medit Link Express. Switch to Medit Link for advanced features and app-based workflows.

#### No Extra Fees. Full Access.

Enjoy complete functionality across both platforms without additional subscriptions.





#### **Simple Process Preferred Type**

For clinics new to digital or focused on fast, simplified workflows.

Prioritizing communication and patient engagement.



#### **Advanced Workflow Type**

For clinics experienced in digital dentistry or operating complex environments.



#### **Hybrid & Expansion Type**

For existing Medit users looking to add mobility without disrupting current workflows.

14

# **Medit Apps & Solutions**

The availability of applications may vary depending on regulatory approvals in each country. Certain apps may not be displayed or offered in specific regions due to local regulatory restrictions.

#### Consultation



#### **Scanning Solution**

#### **Design & Tools**

#### **Medit Smile Design**

Visual Smile Simulation for Patient Engagement

- Show realistic before & after images
- Boost case acceptance with personalized smile design









Before & After Simulation

#### **Medit SmartX**

Intuitive All-On-X Scanning Workflow

- Intuitive, step-by-step workflow for easy full-arch implant scanning
- Ensure scan accuracy with real-time library matching and automatic alignment






#### Medit Orthodontic Suite

Chairside Orthodontic Planning & Simulation

- Combine simulation and treatment planning in one app
- Supports multiple treatment types for doctor-led ortho care





#### **Denture Workflow**

Clinically optimized denture scanning workflow


- Effortless workflows with expanded FOV
- Reline scanning and all-in-one functional impression capture



#### **Impression Scanning**

Fast and Accurate Scanning of Traditional Impressions

- Automatically converts impression scans into 3D models
- Supports various impression types for flexible use



#### **Medit ClinicCAD**

Chairside Prosthetic Design, Made Simple

- Supports crowns, inlays, onlays, bridges, and more
- Complete designs in just a few clicks with step-by-step automation





#### **Medit Model Builder**

Fast, Automated Digital Model Creation

- Quickly generate printable models with auto-design features
- $\, \bullet \,$  Includes tools for articulators, attachments, and advanced customization





#### **Medit Splints**

Easy, Custom Splint & Night Guard Design

- Intuitive UI/UX optimized for clinical workflows
- Versatile tools for a wide range of splint designs and indications




#### **Medit Occlusion Analyzer**

Quantitative Occlusion Analysis for Precise Adjustments

- Visualize pressure distribution and occlusal force mapping
- Reduce subjectivity and support accurate diagnosis with measurable data






#### **Medit Design**

Versatile Multi-Tool for Digital Dentistry



• Expand clinical capabilities with flexible, all-in-one features



16

# Medit i900 Mobility Specifications

Category		Medit <i>i</i> 900 Mobility	
Classifications for Protection Against Electric Shock		Class II, Type BF Applied Parts	
Туре		Wireless	
Handpiece	Dimension	209 x 42 x 42 mm (excluding reusable tip and rechargeable battery)	
	Weight	242 g $\pm$ 5% (excluding reusable tip and rechargeable battery)	
	Wireless Module	Wi-Fi 6E, Bluetooth	
Reusable Tip	Dimension - Full Tip	Large 36.1 x 34.1 x 90.8 mm (W x H x L)	
		Medium 36.1 x 34.1 x 90.4 mm (W x H x L)	
		Small 36.1 x 34.1 x 90.3 mm (W x H x L)	
	Dimension - Tip Head	Large 26.9 x 19.7 mm (W x H)	
		Medium 22.4 x 16.3 mm (W x H)	
		Small 18.36 x 13.1 mm (W x H)	
	Weight	Large 14 g ± 5%	
		Medium 12 g ± 5%	
		Small 11 g ± 5%	
Battery Specification	Working hour	1hr	
	Charging hour	About 2hrs	

- Product availability, regulatory approval status, and release timelines may vary by country. For accurate and up-to-date information, please contactyour local authorized sales representative.
- This document is not affiliated with or endorsed by Apple Inc. iPad and iPadOS are trademarks of Apple Inc., used here for reference. For official information, visit apple.com

#### **About Us**

At Medit, we go beyond creating products to deliver infinite possibilities for our users.

Our commitment to innovation and customer care propels the digital dentistry revolution, offering 3D digital solutions that transform dental clinics and labs.

With our award-winning intraoral scanners, based on our own patented state-of-the-art technology and collaborative software, we are setting new standards in dental technology.

## **Our Medit Community**

With continuous feedback, active discussion, and helpful suggestions, the Medit community is always connected as customers and colleagues actively support and help one another. There's a reason why over 100,000 dental professionals around the globe trust us.

Scan the QR code to join the Medit revolution worldwide and experience it yourself!





