

## Bariatric Steel Wheelchair



This attendant chair was designed to cater for the growing bariatric product market.

### Product features

- Dual, reinforced steel cross braces give extra support and stability
- Generous weight capacity of 200kg (32st)
- Reinforced steel frame provides added support
- Heavy-duty, nylon reinforced upholstery with a back carry pouch for storage
- Removable, reversible desk length armrests facilitate easy side transfers and allow unhindered access to work surfaces
- Height adjustable footrests to suit user support requirements
- 12" rear solid tyres are durable and low maintenance
- Stepper tube allows attendant to safely mount curbs

CODE	WIDTH (OPEN)	WIDTH (CLOSED)	LENGTH	HEIGHT	SEAT WIDTH	SEAT LENGTH	SEAT HEIGHT	WEIGHT	WEIGHT CAPACITY
BTR22BLKAU	67cm (27")	29cm (12")	97cm (38")	90cm (36")	55cm (22")	46cm (18")	52cm (21")	22kg (49lb)	200kg (32st)

## Heavy Duty Rollator

This bariatric rollator is suitable for users up to a higher weight capacity of 180kg (28st), and has been developed to provide the larger person with increased sense of security and safety.

### Product features

- Easy to operate walking aid developed for the larger person to increase their sense of security and safety
- Suitable for bariatric use with a higher weight capacity 180kg (28st)
- Anatomically designed handles allow comfortable grip
- Arthritic friendly looped cable brakes activate effortlessly
- Suitable for both indoor and outdoor use
- Padded seat flips up easily for access to storage basket or for folding
- Padded backrest offers comfortable support
- Foldable for easy storage and transportation
- Handle height can be adjusted between 88cm - 100cm (34.5" - 39") to suit the user
- Fitted with solid 20cm (8") puncture proof tyres for easy manoeuvrability and durability
- Includes a basket for storing essential items
- Available in a Metallic Blue finish



CODE	WIDTH	DEPTH	HEIGHT MIN	HEIGHT MAX	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT MIN	SEAT HEIGHT MAX	WEIGHT	WEIGHT CAPACITY
10215BL-1AU	76cm (30")	69cm (27")	88cm (34.5")	100cm (39")	46cm (18")	33cm (13")	55cm (22")	55cm (22")	12kg (21lb)	180kg (28st)

# Sentra EC Wheelchair *Bariatric*

The Sentra EC Wheelchair has been designed for the larger person to ensure their mobility experience is as comfortable as possible.



- Sturdy, carbon steel frame with durable chrome for an attractive, chip proof, maintainable finish
- Dual axle position allows transition of seat height to hemi level
- Double embossed vinyl upholstery provides comfortable support and is durable, lightweight, attractive and easy to clean
- Simple folding mechanism facilitates storage and transportation
- Footrests can be height adjusted tool free to suit user's support requirements, can be swung out of the way to facilitate transfers or removed for storage/transportation
- Padded desk style armrests provide additional comfort and allow unhindered access to work surfaces
- Armrests are easily removable for sideways transfers, transportation and storage
- Solid tyres, mounted on composite mag rear wheels provide durability and low maintenance
- 20cm (8") solid front castors also ensure longevity and minimum maintenance and are height adjustable in 3 positions
- Maintenance-free, sealed, precision bearings
- Plastic hand rims enable easy manoeuvring for the independent user
- Attendant and user accessible parking brakes on rear wheels provide additional safety
- Elevating footrests are available for users who need to raise their legs, keep them straight or who require a more comfortable foot support
- Available in 3 different seat-widths to suit varying user sizes
- Extra wide seat-widths with reinforced framework are also available – 66cm (26"), 71cm (28"), 76cm (30") with a higher weight capacity of 318kg (50st) Please enquire for details

## Generous weight capacities

CODE	DESCRIPTION
STD20ECDDA-SF	Sentra EC Wheelchair 51cm (20") Self Propel
STD22ECDDA-SF	Sentra EC Wheelchair 56cm (22") Self Propel
STD24ECDDA-SF	Sentra EC Wheelchair 61cm (24") Self Propel
STD26DFA-SF	Sentra EC Heavy Duty Wheelchair 66cm (26") Self Propel
STD28DFA-SF	Sentra EC Heavy Duty Wheelchair 71cm (28") Self Propel
STD30DFA-SF	Sentra EC Heavy Duty Wheelchair 76cm (30") Self Propel



CODE	WIDTH (OPEN)	WIDTH (CLOSED)	LENGTH	HEIGHT	SEAT WIDTH	SEAT LENGTH	SEAT HEIGHT	WEIGHT	WEIGHT CAPACITY
STD20ECDDA-SF	75cm (30")	32cm (13")	107cm (42")	91cm (36")	51cm (20")	41cm (16")	50cm (20")	22kg (49lbs)	200kg (32st)
STD22ECDDA-SF	80cm (32")	34cm (14")	107cm (42")	91cm (36")	56cm (22")	46cm (18")	50cm (20")	26kg (57lb)	200kg (32st)
STD24ECDDA-SF	85cm (34")	37cm (15")	107cm (42")	91cm (36")	61cm (24")	46cm (18")	50cm (20")	26kg (57lb)	200kg (32st)
STD26DFA-SF	89cm (35")	37cm (15")	119cm (47")	97cm (38")	66cm (26")	50cm (20")	43cm-48cm (17"-19")	40kg (88lb)	318kg (50st)
STD28DFA-SF	94cm (37")	37cm (15")	119cm (47")	97cm (38")	71cm (28")	50cm (20")	43cm-48cm (17"-19")	40kg (88lb)	318kg (50st)
STD30DFA-SF	99cm (39")	37cm (15")	119cm (47")	97cm (38")	76cm (30")	50cm (20")	43cm-48cm (17"-19")	40kg (88lb)	318kg (50st)