




Power supply unit/device compatibility list

English

This supplement specifies the power supply units (PSU) compatible with ResMed devices. For full details regarding the correct use of these devices, please see the relevant device User Guides.

Notes:

- Not all products are available in all regions.
- The PSU part number is located in the label at the back of the unit.


Device	Power supply unit						
	Model number	370001	370006	370006	369102	R270-7198	370005
	Part number	-	R370-7555	R370-7565	-	-	-
							
		Original PSU for Air10™ and Lumis devices.	Compact PSU for Air10 and Lumis S/ST devices	Compact PSU for Air10 and Lumis devices	Packed with the S9™ series.	Original PSU for Astral.	Compact PSU for Astral series.
AirSense™ 10 / AirStart™ 10		✓	✓	✓	-	-	-
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST		✓	✓	✓	-	-	-
AirCurve™ 10 ASV / AirCurve™ 10 ASVAuto		✓	✓	✓	-	-	-
AirCurve™ 10 CS Pacewave™		✓	✓	✓	-	-	-
Lumis™ S / Lumis™ ST		✓	✓	✓	-	-	-
AirCurve™ CS-A		✓	-*	✓	-	-	-
AirCurve™ ST-A		✓	-*	✓	-	-	-
Lumis™ ST-A		✓	-*	✓	-	-	-
S9™		-	-	-	✓	-	-
ResMed Power Station II		✓	-*	-*	✓	-	-
ResMed Power Station II-V2		✓	✓	✓	-*	-	-
Stellar™		-	-	-	-	-	-
Astral™		-	-	-	-	✓	✓

*These devices have not been validated for use with this PSU. The PSU will deliver the same power to the unit and meet the power requirements of the device, however the Electromagnetic Compatibility has not been verified.

⚠ WARNING

Only use compatible ResMed power supplies with ResMed devices.

Sample label:



Model number

Part number