

Anti-Draft Chamber 153



More information on the website radwag.com/en/info,w1,EL5



Overview


The anti-draft chamber is intended for balances featuring a weighing pan 128×128 mm. It is an optional equipment for balances operating in unfavourable ambient conditions (air movement air drafts). It is installed as replacement of a standard pan draft shield. The anti-draft chamber features slidable top and side doors ensuring unbounded access to the weighing pan.

Compatible with


PS X2 Precision Balance PS R1 Precision Balance PS R2 Precision Balance


Device dimensions


WLC X2 Precision Balance PS X2 Balances for PGC














Anti-draft chamber for X2 and R2 precision balances







Anti-draft chamber for balances

with 128x128 mm weighing pan







Anti-draft chamber for APP-00-SK mass comparators






Anti-draft chamber for WAY-KO-SK mass comparators








80,5

Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for XA 4Y / XA 4Y PLUS balances