

Focus on Liquid Handling





One company, *four internationally-recognised brands*

Introduction

Stuart offers a comprehensive range of liquid handling products. The Bibbyjet Pro is a powerful and light weight rechargeable pipette controller. The Bibbypet range of bottle top dispensers offer an economical solution to measured bottle dispensing, with the Pressmatic range available for more precise dispensing. The range is completed with the digital burette, making titrations quick and easy.



>>

Bibbyjet Pro Pipette Controller

- Dual-controls provide simple operation for filling and delivery
- One full working day of pipetting without recharging
- Specially designed silicone adapter provides a tight fit for all glass or plastic pipettes from 0.1 to 100ml
- Built-in non-return value combined with a membrane filter effectively protects against penetration of liquids

Digital Burette

- Accuracy of <0.2%, Precision of < 0.1%, Resolution of 0.01ml
- Lifetime battery for over 60,000 titrations with 'auto-off' function
- 'Zero' reset button allows for quick and simple completion of multiple titrations
- Instant and accurate large LCD read-out provides constant display of volume dispensed

Pressmatic HF Bottle Top Dispenser >>

- High precision dispensing for highly aggressive liquids such as hydrofluoric acid
- Piston, barrel , valve block and tubes are manufactured from chemically inert ceramic and PTFE
- Innovative ergonomic design that maximises safety and ease of use
- Recirculation valve delivers the reagent used for priming back into the bottle
- Available in 2.0 to 10.0ml capacity

Bibbypet Bottle Top Dispenser >>

- Economical solution for the dispensing of non-aggressive liquids
- Accurate and reliable dispensing of many common laboratory liquids
- Low-cost replaceable piston/cylinder and valves
- Available in either:
 - 1 to 10ml
- 2.5 to 25ml



Pressmatic

10-60



< Pressmatic PP Bottle Top Dispense

- Accuracy of <0.7% and Repeatability of <0.1%
- Accurate dispensing of many common laboratory liquids
- Simple click stop ring for fast and accurate volume setting
- Clear viewing window to check priming
- Available in either:
 - 2.0 to 10.0ml
- 10.0 to 60.0ml

Pressmatic R Bottle Top Dispenser >>

- High precision dispensing ideal for most solutions, including concentrated acid and alkali solutions
- Innovative ergonomic design that maximises safety and ease of use
- Recirculation valve delivers the reagent used for priming back into the bottle
- Chemically resistant and non-sticking ceramic piston
- Available in either:
 - 0.2 to 1.0ml
 - 0.4 to 2.0ml
 - 1.0 to 5.0ml
 - 2.0 to 10.0ml
 - 5.0 to 30.0ml
 - 10.0 to 60.0ml

You Tube

Bibby Scientific videos now on YouTube! www.youtube.com/bibbyscientific

Pressmatic	: Bottleto	p Dispen	sers						Ŧ
Model	PD10P	PD60P	PD1R	PR2R	PD5R	PD10R	PD30R	Burnatik D-60w V⊇H PD60R	PD10R/HF
Range	Pressmatic	Pressmatic	Pressmatic	Pressmatic	Pressmatic	Pressmatic	Pressmatic	Pressmatic	Pressmatic
Volume, ml	2.0 to 10.0	10.0 to 60.0	0.2 to 1.0	0.4 to 2.0	1.0 to 5.0	2.0 to 10.0	5.0 to 30.0	10.0 to 60.0	2.0 to 10.0
Subdivision, ml	0.25	1.00	0.05	0.05	0.10	0.25	0.50	1.00	0.25
Accuracy, %	<0.7	<0.7	<0.7	<0.7	<0.7	<0.7	<0.7	<0.7	<0.7
Repeatability, %	<0.1	<0.1	<0.2	<0.1	<0.1	<0.1	<0.1	<0.1	<0.1

Bibbypet		
Model	BD10	BD25
Range	Bibbypet	Bibbypet
Volume, ml	1 to 10	2.5 to 25
Subdivision, ml	0.2	0.5
Accuracy, %	1.2	1.2

Bibbypet dispensing cartridge

Model	Capacity	Materials	Pack Qty
BD1010	10ml	Piston PE, cylinder PP	3
BD1012	25ml	Piston PE, cylinder PP	3

Bibbypet replacement parts

Model	Description	Pack Qty
BD1014	Discharge tube	1
BD1016	Filling tube, 250mm	2

Polypropylene adapters

Model	Description	Thread
PD4532PP	Polypropylene bottle top adapter	32mm
PD4538PP	Polypropylene bottle top adapter	38mm
PD4540PP	Polypropylene bottle top adapter	40mm
PD4532PT	Chemically resistant ECTFE bottle top adapter	32mm
PD4538PT	Chemically resistant ECTFE bottle top adapter	38mm
PD4540PT	Chemically resistant ECTFE bottle top adapter	40mm

BibbyJet Pipette Controller

	•
PC2000	Bibbyjet Pro with charger for
	UK and Ireland
PC2001	Bibbyjet Pro with charger for
	Europe
PC1011	Silicone Adapter
PC1012	0.2µm filter, sterile, pack of one
PC1013	0.2µm filter, sterile, pack of ten
PC1014	3µm filter, non-sterile, pack of one
PC1015	3µm filter, non-sterile, pack of ten
PC1016	Wall support

Bibby Scientific

Bibbyjet Pro

Digital Burette, BT30, BT50

Model Description BT30 30ml Digital Burette with 33, 38 and 45mm adapters BT50 50ml Digital Burette with 33, 38 and 45mm adapters

Accessories

Model	Description
BTA33	Bottle adapter, 33mm thread
BTA38	Bottle adapter, 38mm thread
BTA40	Bottle adapter, 40mm thread
BTA45	Bottle adapter, 45mm thread

Digital Burette Reservoirs

Description	Capacity
Reservoir with side filling port chamber	1400ml
Reservoir with side filling port chamber	2500ml
	Reservoir with side filling port chamber



Bibby Scientific - UK Beacon Road, Stone, Staffordshire, ST15 0SA, United Kingdom T: +44 (0)1785 812121 F: +44 (0)1785 813748 E: sales@bibby-scientific.com

Bibby Scientific - France Bâtiment le Deltaparc Silic Paris Nord 2, 7 rue du Canal, BP 55437 Villepinte 95944 ROISSY Charles de Gaulle, France T: +33 (0) 1 48 63 78 00 F: +33 (0) 1 48 63 78 01 E: ventes@bibby-scientific.fr

Bibby Scientific - US 3 Terri Lane, Suite 10, Burlington, NJ 08016, USA T: +1 609 589 2560 F: +1 609 589 2571 E: labproducts@techneusa.com

Bibby Scientific - Singapore Prudential Tower, Level 26 30 Cecil Street, Singapore, 049712 T: + 65 6631 2976 F: +44 (0) 1785 810405 E: info@bibby-scientific.com

Bibby Scientific - Middle East PO Box 27842 Engomi 2433 Nicosia, Cyprus T: + 357 22 660 423 F: + 357 22 660 424 E: sales@bibby-scientific.me