

Digital Cel-Gro Tissue Culture Rotator

Brushless no-maintenance motor
provides a gentle rotating motion for
optimum results



The Thermo Scientific™ Digital Cel-Gro Tissue Culture Rotator precisely mixes test tube contents. Achieve optimum culture results for tissue and virus growth, hormone production, and sample testing of viruses and cytotoxicity assays.

Key Features

- Digital speed and time display
- Adjustable angle supports mixing requirements for a diverse range of applications
- Variable-speed permits a wide range of motion
- Brushless, no maintenance motor provides smooth operation and worry-free use
- Drum selection accommodates various capacities
- Speed range: 5 to 70rpm
- Speed accuracy: ±1rpm
- Adjustable angle: 95°- 180°
- Timer range: up to 166 hours
- Maximum load: 5kg (11 lbs.)

Drum sold separately (required for operation)
Warranty: 2 years | Certifications: cCSAus, CE



Digital Cel-Gro Tissue Culture Rotator assembled with Thermo Scientific Rotator Drum and Nunc Test Tubes. Nunc Test Tubes sold separately.

Product	Electrical Requirements	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digital Cel-Gro Tissue Culture Rotator	100 to 240V, 50/60Hz, 0.2A	38.3 x 26 x 52.2cm (15.1 x 10.2 x 20.6 in.)	9.5kg (20.9 lbs.)	US	88882015
				EU/UK/CHN	88882016

Accessories

Product	Overall Dimensions	Quantity	Catalog Number
Rotator Drum, 50 x 1.5/2.0mL	Ø29 x 2.4cm (11 x 0.9 in.)	Each	88882149
Rotator Drum, 30 x 15mL	Ø29 x 3.1cm (11 x 1 in.)	Each	88882150
Rotator Drum, 15 x 50mL	Ø29 x 4.8cm (11 x 1.9 in.)	Each	88882151
Rotator Drum, 10 x 10/15mL + 20 x 5/7mL	Ø29 x 3.1cm (11 x 1 in.)	Each	88882152
Rotator Drum, 64 x Ø14mm	Ø31 x 7.5cm (12 x 3 in.)	Each	88882153
Rotator Drum, 64 x Ø18.5mm 15-20mL	Ø31 x 11.5cm (12 x 4.5 in.)	Each	88882154
Rotator Drum, 142 x Ø17mm	Ø38 x 7.8cm (15 x 3 in.)	Each	88882155
Rotator Drum, 76 x Ø26mm	Ø38 x 7.3cm (15 x 2.9 in.)	Each	88882156
Rotator Drum, 60 x Ø30mm	Ø38 x 7.8cm (15 x 3 in.)	Each	88882157

Product	Overall Dimensions	Quantity	Catalog Number
Power Adaptor, US plug (125VAC, 10A)	1.8m	Each	88870126
Power Adaptor, AUS/CHN plug (250VAC, 10A)	1.8m	Each	88870127
Power Adaptor, EU plug (250VAC, 16A)	1.8m	Each	88870128
Power Adaptor, UK plug (250VAC, 13A)	1.8m	Each	88870129



Digital Cel-Gro Tissue Culture Rotator Drums, 38cm diameter.