Thermo Scientific Owl A6 Horizontal Electrophoresis System

Ideal for high-throughput laboratories



- System screens 25 to 500 samples on a single agarose gel, producing clear, tight banding patterns with no "smiling"
- Gel tray's unique comb-slot design allows use of up to 10 evenly spaced combs
- With special combs and a multichannel pipetter, samples are loaded directly from a 96-well plate, 8 or 12 at a time
- Incorporates unique end-gate casting design gel tray is fitted with gasketed end gates which provide a leakproof seal without tape
- Sample loading with multichannel pipetter speeds up transfer
- Bubble level ensures casting of flat, uniformly thick gels

Includes: Buffer chamber with three-point leveling base, SuperSafe* lid with attached power supply leads, UVT gel tray with gasketed end gates, four combs (two 25-well and two 50-well, 1.5mm thick)

Warranty¥: Three years

Cat. No.	Horizontal System Description	Gel Size (L × W)	Footprint (L × W × H)	Running Buffer Volume	Sample Range
A6	Wide	25 × 23cm	35 × 38 × 11cm (13.75× 14.94 × 4.38in.)	3.0L	25 to 500

Thermo Scientific Owl A6 Wide Gel Electrophoresis System Replacement Parts

- Power supply leads include a retractable sheath end for convenience and safety
- Power supply leads are manufactured with strain relief cord to prevent breakage
- UVT gel tray incorporates gasketed end gates which provide a leakproof seal without tape

Cat. No.	Description
PSL-5	Power supply leads, 3 ft. (Includes one black and one red lead)
A6-UVT	UVT gel tray with gasketed end gates
A3-GK-1	Replacement gaskets (1 pair)
A3-1-EG Replacement end gates (1 pair)	

24



Thermo Scientific Owl Leveling Platform

- Three-point leveling platform ensures flat casting and running surface
- Includes bubble level

Cat. No.	Туре
BBL-1	Bubble Level



Thermo Scientific Owl Casting Dams

- Seal agarose on contact, eliminating the need for tape
- Do not attach to gel tray, enabling gels of any length to be poured
- Easily handled and maneuvered using knob or finger depressions
- High-quality aluminum dam blocks off a portion of UVT gel tray, allowing shorter gels to be cast and run in one device

Cat. No.	Description
DAM-23	Chameleon casting dam; 23cmW

Thermo Scientific Owl A6 Horizontal Electrophoresis System Combs

Cat. No.†	Description	Thickness / Width of Teeth	RECOMMENDED LOADING VOLUMES Calculated at 75% of total well volume			
			0.25cm gel thickness	0.5cm gel thickness	0.75cm gel thickness	1.0cm gel thickness
D3-MTC‡	Microwell 25 (1X)-tooth comb	1.0mm / 7.5mm	6µL	20µL	34µL	48µL
D3-MTD‡	Microwell 25 (1X)-tooth comb	1.5mm / 7.5mm	8µL	30µL	51µL	72µL
D3-MT2C‡	Microwell 50 (2X)-tooth comb	1.0mm / 3mm	2μL	8µL	14µL	19µL
D3-MT2D‡	Microwell 50 (2X)-tooth comb	1.5mm / 3mm	3µL	12µL	20μL	29μL
D3-WALL	Wall comb, 1 tooth	1.5mm / 230mm	_	_	_	_

[†] Customized combs are available (Order SKU # XCM). Please contact your Customer Service Representative for further details.

^{‡ 8} and 12-channel pipette format (1X) and (2X) combs.