

Whitley A55 Anaerobic Workstation - Porthole Options



Instant Access Porthole System / Instant Access Porthole System

This version of the Whitley A55 Anaerobic Workstation is equipped with 4x Instant Access Porthole Systems.

This version is ideal for users wishing to enter the workstation frequently for brief working periods.

**Option 1 can be ordered with factory modifications to accept conventional oval, sleeved ports as part of a replacement removable front. For this version please order:
230V: Order A07710 / 110V: Order A07711**

230V: Order A07000 / 110V: Order A07001



Instant Access Porthole System / Conventional Porthole System

This version of the Whitley A55 Anaerobic Workstation is equipped with 2x Conventional Porthole Systems with sleeves and cuffs on the right-hand side, on the removable front and 2x Instant Access Ports on the left-hand side.

230V: Order A07720 / 110V: Order A07721



Conventional Porthole System / Conventional Porthole System

This version of the Whitley A55 Anaerobic Workstation is equipped with 4x Conventional, oval ports.

This version is ideal for users wishing to enter the workstation less frequently for longer working periods.

230V: Order A07730 / 110V: Order A07731



Conventional Porthole System / Instant Access Porthole System

This version of the Whitley A55 Anaerobic Workstation is equipped with 2x Conventional Porthole Systems with sleeves and cuffs on the left handside and 2x Instant Access Ports on the right-hand side.

230V: Order A07740 / 110V: Order A07741

Note: All options of the Whitley A55 have a removable front, on the right-hand side, covering an opening which measures 800mm (31.5") wide x 398mm (15.5") high