

FACT SHEET

WHITLEY A45 WORKSTATION

The A45 Anaerobic Workstation has <u>three</u> gloveless/sleeveless ports for convenient user access to the entire incubation and working areas.



Features and Benefits

- Accommodates over 750 x 90mm Petri dishes whilst still providing a generous working area.
- Features a revolutionary new porthole system that allows entry to the chamber in seconds and needs neither gloves nor sleeves, which means that sample manipulation within a secure anaerobic environment is now as easy as working aerobically on the laboratory bench.
- Optional removable front (with instant access or traditional sleeved ports) to enable the transfer of larger batches of samples or pieces of equipment for use inside the chamber.
- Built-in rapid airlock allows the transfer of plates in just 60 seconds.
- Colour, touch-screen control panel for ease of use and for visual display of parameters such as temperature, humidity, and airlock cycle status.
- Available with fully integrated Anaerobic Conditions Monitor with data download facility.
- Automatic commissioning cycle to save you time and resources.
- Optional fully comprehensive maintenance and breakdown contracts are available to prolong the life of your investment. (Note: different arrangements apply in different countries).
- DWS provides scientific back-up from our fully operational microbiology laboratory (UK only).

Facts and Figures

Weight: 132kg Dimensions: 720 x 1660 x 710mm (D x L x H) Gas Requirements: 1 x cylinder of Anaerobic Gas Mix and 1 x cylinder of Nitrogen Electricity Requirements: 230 \pm 10% V AC or 115 \pm 10% V AC – single phase