



***Mycoplasma bovis* Antibody Test Kit (ELISA)** **Bovichek® *M. bovis* – 1 or 2 plates**

Recently IMPROVED *M. bovis* Antibody Test Kit offers the possibility to test milk samples in addition to serum samples

Materials provided with the kit (2 plates)

Components	Quantity
• 12 strips of 8 wells coated with <i>M. bovis</i> Ag	2
• Ready-to-use positive control	6 mL
• Ready-to-use negative control	6 mL
• Ready-to-use conjugate	25 mL
• Concentrated wash solution (10X)*	2 x 125 mL
• Ready-to-use substrate	30 mL
• Ready-to-use stop solution	15 mL

* Crystals may form when wash solution is kept at 2-7°C. This will not affect the efficiency of the product. In order to use this solution, simply bring it to room temperature and the crystals will dissolve.

Materials Required but not Provided

- Purified water
- Adjustable single and multi-channel micropipettes
- Single-use micropipette tips
- ELISA microplate washer
- Test tubes for sample dilution
- ELISA 96-well microplate reader equipped with 450 nm filter
- Containers for dilution of other solutions

The materials provided are sufficient for testing up to 184 samples.

IMPORTANT:

- Requires a permit to release veterinary biologics from CFIA in Canada. This product is not licensed by CFIA and any claims made have not been substantiated by CFIA.
- Requires a Research and Evaluation import permit from the USDA in the US. This product is not licensed by the USDA and any claims made have not been substantiated by the USDA.

Rév. 2017-06-19

For further information visit www.biovet-inc.com or contact customer service:

USA
1-877-8BIOVET
 (824-6838)



Biovet®
order@biovet-inc.com

International & Canada
1-888-8BIOVET
 (824-6838)