



OVUCHECK® Plasma

Measurement of Progesterone in Plasma or Serum

THE PRECISION OF A REFERENCE METHOD AT A MUCH LOWER COST!

The OVUCHECK® PLASMA ELISA kit is an immunoenzymatic test which provides a simple, reliable and precise measurement of progesterone in plasma or serum of cows and bitches.


OVUCHECK® PLASMA is used for oestrus detection and assessment of pregnancy status and luteal function.

Major advantages of OVUCHECK® Plasma

- No more waiting until the next day for results
- Reliable results within 85 minutes with an ELISA reader
- Precision comparable to that of a reference method
- Only requires an ELISA reader
- Much lower cost analyzer
- Much lower cost reagent
- Precise progesterone analyses made more affordable than ever!

Kit components	Quantity
• 12 strips of 8 wells coated with progesterone antibodies	1
• Ready-to-use progesterone standards (Each of: 1.0; 2.5; 5.0; 10.0 ng/mL)	1.0 mL
• Ready-to-use conjugate-(Prog.-AP)	26 mL
• Substrate buffer	25 mL
• Substrate tablets	3 X 40 mg
• Ready-to-use stop solution	20 mL
• Graph paper	1

Each kit contains sufficient reagents for up to 92 tests plus four standards.



Materials required but not Provided:

- Purified water
- Adjustable single- and multi-channel micropipettes
- Single-use micropipette tips
- ELISA 96-well microplate reader equipped with 405 nm filter

For further information visit www.biovet-inc.com or contact customer service:

USA
1-877-8BIOVET
(824-6838)

Biovet®
order@biovet-inc.com

International & Canada
1-888-8BIOVET
(824-6838)