

Composition**P04-48010P GBSS**
w/o: NaHCO₃

	Components	mg/L
Inorganic Salts	Calcium chloride x 2H ₂ O	225.10
	Magnesium chloride x 6H ₂ O	210.00
	Magnesium sulfate anhydrous	34.20
	Potassium chloride	370.00
	Potassium dihydrogen phosphate anhydrous	30.00
	Sodium chloride	7,000.00
	di-Sodium hydrogen phosphate	120.00
Other Components	D(+)-Glucose anhydrous	1,000.00