


## IVF & insemination syringe 1 mL

1. Company profile		Update : November 2016 Issue date : August 2015
1.1	<b>Name :</b> Laboratoire CCD	
1.2	<b>Address :</b> 48, rue des Petites Ecuries 75010 Paris FRANCE	<b>Tel :</b> + 33 (0)1 44 95 14 95 <b>E-mail :</b> contact@ccdlab.com <b>Site internet :</b> www.ccd-lab.com
1.3	<b>Medical Device Vigilance contact :</b> Emilie FEZELOT, Q/A manager	<b>Tel :</b> +33 (0)3 44 26 89 56 <b>E-mail :</b> efezelot@ccdlab.com

2. Product information	
2.1	<b>Common designation :</b> 1mL syringe
2.2	<b>Trade name :</b> IVF & AI Syringe 1mL
2.3	<b>Medical device classification :</b> Class IIa <b>Relevant EU Directive :</b> 93/42/CEE <b>Certified Body identification number :</b> CE 0459 <b>Manufacturer :</b> Laboratoire CCD
2.4	<b>Product features :</b> This syringe is designed specifically for Medically Assisted Procreation procedures. It is composed of three parts: <ul style="list-style-type: none"> <li>▪ A body in polypropylene</li> <li>▪ A plunger in polypropylene</li> <li>▪ A seal in polyisoprene</li> </ul> Volume: 1 mL
2.5	<b>PRODUCT CODE :</b> 132310 <b>Packaging :</b> Unit packaging <b>Ordering unit :</b> Unit <b>Minimum order :</b> box of 25 units <b>Labelling :</b> <div data-bbox="231 1565 884 1881" data-label="Image">  </div>

<b>2.6</b>	<b>Components and accessories :</b> Polypropylene  <b>Active substances :</b> Latex free Phtalates free Contains no component of biological or animal origin
<b>2.7</b>	<b>Indications :</b> <ul style="list-style-type: none"> <li>▪ Artificial Insemination</li> <li>▪ Embryo transfer and in vitro fertilization</li> </ul>

### 3. Sterilization method

	<b>Sterile medical device :</b> Yes <b>Sterilization method :</b> Irradiation
--	--

### 4. Storage and shelf life

	<b>Current conservation and storage conditions :</b> Store in a dry place, away from light <b>Specific precautions :</b> Do not use if unit packing is damaged <b>Product shelf life :</b> 5 years
--	---

### 5. Operational safety

<b>5.1</b>	<b>Technical safety :</b> Single use medical device. Do not re-sterilize. Do not re-use.
<b>5.2</b>	<b>Biological safety :</b> N/A

### 6. Instructions for use

<b>6.1</b>	<b>Operating instructions :</b> The 1 mL syringes for IVF and AI are designed to be connected to insemination catheters and embryo transfer catheters.
<b>6.2</b>	<b>Indications :</b> <ul style="list-style-type: none"> <li>▪ Artificial Insemination</li> <li>▪ Embryo transfer and in vitro fertilization</li> </ul>
<b>6.3</b>	<b>Precautions for use :</b> <ul style="list-style-type: none"> <li>• Check that the packaging is not damaged</li> <li>• Before use, check the integrity of the medical device</li> <li>• Prior to use, rinse the device with appropriate medium according to state of the art practices</li> </ul>
<b>6.4</b>	<b>Contra-indications :</b> <ul style="list-style-type: none"> <li>• None</li> </ul>

**7. Additional information on the product**

N/A

**8. Appendices**

N/A

**9. Pictures**
