

# Tabletop Autoclave

## HMT 232N



The HMT 232 N is an easy-to-use, semi-automatic autoclave with high performance, reliability and safety.



For use in non medical laboratories for sterilizing liquids, objects, etc. The chamber consists of easy-to-clean stainless steel.

Different programs to choose from.

The sterilization pressure is adjustable from 0.9 to 2.1 bar.

For the sterilization following settings are available:

4 minutes / 15 min / 30 min.

equipped with digital sterilisation cyler recorder (SD & USB)

Automatic control by the end of the sterilization process, the entire sterilization process is microprocessor controlled and monitored

### SAFETY:

- It features a low water protection, safety and venting valve, pressure switch and overtemperature protection.
- cover-impact backup
- Each element is equipped with an independent thermal overload switch.
- CE certified

MODEL	HMT-232N
Chamber dimensions (mm):	230mm Ø X 410mm(D)
Chamber volume	ca. 16Liter
External dlemension:	501(W) X406(H) X 537(D)mm
Power supply:	230V / 1,4KW, (max.) 6,5A
Working temperature:	115-134°C; 2.1bar
Max. temp.:	132°C
Safety:	Overheat protection Pressure safety valve Over pressure lock Micro switch door lock
Accessories(optional)	Sterilisation boxbox Tray set



Kellerstraße1  
D-84577 Tüßling  
Tel.: +49(0)8633/50 52 0-0  
Fax.:+49(0)8633/50 52 0 - 99  
[info@hmc-europe.com](mailto:info@hmc-europe.com)  
[www.hmc-europe.com](http://www.hmc-europe.com)