



DISTRIWEL 440[®]

Media Pourer Stacker





DISTRIWEL 440®

Increase your productivity
quality
traceability

PERFORMANCES



440 plates
poured
in 30 mins

DISTRIWEL 440®

Walkaway device!

Up to 440 dishes
filled in 30 mins



0°



Ensure a flat poured media
(bubble level integrated)

7" touch screen LCD
Intuitive display

Traceability
Print up to 25 characters
(2 lines) on the Petri dish
(optional laser printer)



The Built-in Peltier cooling
effect reduces condensation
and solidification time

AEROBIOCONTAMINATION PREVENTION

Contamination prevention
(UV lamp) while pouring
+ Plexiglas protection cover

In-line blood
addition option

Multicolour
visual alarms

Your culture media is accurately
poured by a Masterflex™
peristaltic pump



SPECIFICATIONS

DISTRIWEL 440®



DISTRIWEL FEATURES

Dishes/hour for 90mm	Up to 750
Capacity of the carousel	440 dishes
Fast plate cooling	Built-in Peltier effect
Contamination prevention	Built-in UV lamp
Pour plates function	✓
Pouring accuracy	±1%
Pouring range	1 – 99 ml
Stack plates function	✓
Peristaltic pump	Integrated Masterflex pump
Stored programs	40
Control panel	LCD Touch screen
Country of origin	100% France
Traceability (optional)	Laser printer (LINX CJ400) → Printed Petri dishes
Footswitch (for manual dispensing)	✓ (optional)

DIMENSIONS & ELECTRICAL SUPPLY

W x H x D	808 x 940 x 689 mm
Electrical supply	100-250v 50-60 Hz 10 A
Power supply	150 W

PACK YOUR
READY
TO USE DISHES !



ORDER IT TOGETHER WITH DISTRIWEL 440 : Ref. MINIMA



ALLIANCE BIO EXPERTISE



Empowered from our people's skills and experiences, we, ALLIANCE BIO EXPERTISE, together with our customers and partners decided to manufacture Microbiology automated equipment around essential values:
Customer Focus, Service, Loyalty, Innovations

**Ask your local partner
for a personal demonstration!**

ZI de courbouton , Le Tremplin – 35480 GUIPRY - FRANCE
E-mail : contact@abioexpertise.com
Phone : +33(0)2.40.51.79.53
alliance-bio-expertise.com