
Acetylthiocholin-jodid reinst Acetylthiocholine-iodide research grade	Kat.Nr./Cat.No. :	10570
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_7H_{16}N^+OS \cdot I^-$
Molekulargewicht molecular weight		289.2
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	HPLC	min. 98.0
Löslichkeit solubility	0.1 M in H ₂ O	klar und farblos clear and colorless
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.