

Cycloheximid krist. rein Cycloheximide cryst. pure	Kat.Nr./Cat.No. :	10700
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₅ H ₂₃ NO ₄
Molekulargewicht molecular weight		281.4
Aussehen appearance		weißes bis graugelbes Pulver white to grayish yellow powder
Biologische Aktivität (U/mg) biological activity		> 900
Gehalt (%) assay	HPLC	min. 90.0
Schmelzpunkt (°C) melting point		98 – 112
Identität identity	IR	entspricht corresponds
Ursprung origin	nicht tierisch, isoliert aus Streptomyces narases non animal, isolated from Streptomyces narases	
Lagerung (°C) storage		+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--