
Ethylendiamin-tetraessigsäure p.A. EDTA, analytical grade	Kat.Nr./Cat.No. :	11278
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₀ H ₁₆ N ₂ O ₈
Molekulargewicht molecular weight		292.3
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration	min. 99.0
Trocknungsverlust (%) loss on drying	3 h, 105°C	max. 0.5
Identität identity	IR	entspricht corresponds
A 1cm/ 280 nm	5% in 1M NaOH	< 0.25
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.