

---

**Albumin human** **Cat.No. : 11877**  
lyophilized

---

| Parameter                        | Method                   | Specification                  |
|----------------------------------|--------------------------|--------------------------------|
| <b>Appearance</b>                |                          | pale yellow to greenish flakes |
| <b>Albumin purity (%)</b>        |                          | > 97                           |
| <b>Total Proteins (%)</b>        |                          | > 92                           |
| <b>Moisture (%)</b>              |                          | < 5                            |
| <b>pH</b>                        | 10 % in H <sub>2</sub> O | 6.5 – 7.4                      |
| <b>Mesophilic germs (CFU/ g)</b> |                          | < 100                          |
| <b>Storage (°C)</b>              |                          | +2 to +8                       |

The source material was tested and found negative for:

- Antibodies: anti-HIV 1 + 2, anti-HCV
- HBs antigen
- HIV 1 and HCV RNA
- Syphilitic serology

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.