

Borsäure Boric acid electrophoresis grade	Kat.Nr./Cat.No. :	15166
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		H ₃ BO ₃
Molekulargewicht molecular weight		61.83
Aussehen appearance		weißes, kristallines Pulver white, crystalline powder
Gehalt (%) assay	Titration	99.0 – 100.5
pH	Ph.Eur.	3.8 – 4.8
Sulfat (ppm) sulfate	Ph.Eur.	max. 450
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 15
Identität Identity	IR	entspricht corresponds
Anwendung application	DNA-Elektrophorese DNA-electrophoresis	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.