

Caesiumchlorid p.A. Cesium chloride analytical grade	Kat.Nr./Cat.No. :	15554
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		CsCl
Molekulargewicht molecular weight		168.36
Aussehen appearance		weißes krist. Pulver white cryst. powder
Gehalt (%) assay		min. 99.9
A 1cm/ 260nm	50% in H ₂ O	max. 0.02
A 1cm/ 280nm	50% in H ₂ O	max. 0.02
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 1
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.