

<b>Cholsäure-Na-salz, p.A.</b> Cholic acid-Na-salt, analytical grade	<b>Kat.Nr./Cat.No. :</b>	<b>17126</b>
---	--------------------------	--------------

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>24</sub> H <sub>39</sub> O <sub>5</sub> ·Na
<b>Molekulargewicht</b> molecular weight		430.5
<b>Aussehen</b> appearance		weißes Pulver white powder
<b>Gehalt (%)</b> assay	HPLC	min. 98.0
<b>[α]<sub>D</sub>/ 20 °C</b>	10% in H <sub>2</sub> O	+26° bis/ to +28°
<b>Trocknungsverlust (%)</b> loss on drying		max. 5.0
<b>pH</b>	2% in H <sub>2</sub> O	7.5 – 8.5
<b>Schwermetalle (ppm)</b> heavy metals		< 20
<b>Lagerung (°C)</b> storage		+15 bis / to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.