

<b>4',6-Diamidino-2-phenylindol-2HCl p.A.</b>	<b>Kat.Nr./Cat.No. :</b>	<b>18860</b>
4',6-Diamidino-2-phenylindole-2HCl analytical grade		

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>16</sub> H <sub>15</sub> N <sub>5</sub> ·2HCl
<b>Molekulargewicht</b> molecular weight		350.25
<b>Aussehen</b> appearance		gelbes Pulver yellow powder
<b>Reinheit (%)</b> purity		min. 98
<b>Löslichkeit</b> solubility	2.5 % in H <sub>2</sub> O	klare, gelbe Lösung clear, yellow solution
<b>λ<sub>max</sub> (nm)</b>	0.001 % in H <sub>2</sub> O	340 – 348
<b>A 1 cm/ λ<sub>max</sub></b>	0.001 % in H <sub>2</sub> O	min. 0.7
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.