

3,3-Diaminobenzidin·4HCl·xH₂O reinst 3,3-Diaminobenzidine·4HCl·xH ₂ O research grade	Kat.Nr./Cat.No. :	18865
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₂ H ₁₄ N ₄ ·4HCl·xH ₂ O
Molekulargewicht molecular weight		360.1 (wasserfrei/anhydrous)
Aussehen appearance		zartgraues Pulver * light grey powder
Löslichkeit solubility	50 mg/ml H ₂ O	klare dunkelrotbraune Lösung clear dark reddish brown solution
Gehalt (%) assay	Titration	min. 96.0
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

* Das Aussehen dieses Produktes kann produktionsbedingt von hellrosa bis grau-violett variieren.
The appearance of this product may vary due to processing from light rose-colored to grey-violet.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.