

Dodecylsulfat·Na salz in Pellets, reinst Dodecylsulfate·Na salt in pellets, research grade	Kat.Nr./Cat.No. : 20765
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₂ H ₂₅ O ₄ SNa
Molekulargewicht molecular weight		288.38
Aussehen appearance		weiße bis fast weiße Pellets white to off-white pellets
Gehalt (%) assay	Titration	min. 99.5
Löslichkeit solubility	10% in H ₂ O	klar, max. 20 Hazen clear, max. 20 Hazen
Trübung (NTU) turbidity	10% in H ₂ O 20 – 25 °C	max. 2.5
A 1 cm/ 260 nm	10% in H ₂ O	max. 0.2
A 1 cm/ 280 nm	10% in H ₂ O	max. 0.1
Sulfatasche (%) sulfated ash		21.5 – 27.5
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.