
Glutaraldehyd 25% in Wasser Standardqualität Glutaraldehyde 25% in water standard grade	Kat.Nr./Cat.No. :	23115
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₅ H ₈ O ₂
Molekulargewicht molecular weight		100.13
Aussehen appearance		klare, farblose bis gelbliche Lösung clear, colorless to yellowish solution
Gehalt (%) assay		23.0 – 27.0
Brechungsindex refractive index	25°C	1.3690 – 1.3755
pH	25% in H ₂ O	2.5 – 3.0
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.