

<b>Glycerin-2-phosphat Dinatriumsalz Pentahydrat</b> rein 2-Glycerolphosphate-Na <sub>2</sub> -salt Pentahydrate pure	<b>Kat.Nr./Cat.No. :</b>	<b>23330</b>
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		C <sub>3</sub> H <sub>7</sub> O <sub>6</sub> P·Na <sub>2</sub> ·5H <sub>2</sub> O
<b>Molekulargewicht</b> molecular weight		306.1
<b>Aussehen</b> appearance		weiße Kristalle white crystals
<b>Gehalt (%)</b> assay	α + β-Glycero-PO <sub>4</sub>	min. 99.0
<b>Wasser (%)</b> water	KF	26 – 33
<b>Identität</b> identity	IR	entspricht corresponds
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
We do not guarantee that the product can be used for a special application.