
3-Methacryloxypropyltrimethoxysilan 3-Methacryloxypropyltrimethoxysilane (Bind-Silane)	Kat.Nr./Cat.No. :	28739
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₀ H ₂₀ O ₅ Si
Molekulargewicht molecular weight		248.4
Aussehen appearance		klare farblose Flüssigkeit clear colorless liquid
Gehalt (%) assay	GC	98.0 – 100.0
Brechungsindex refractive index	20°C	1.430 – 1.433
Farbzahl color	APHA	max. 30
Viskosität (mm²/s) viscosity		2.3 – 2.5
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.