

L-Milchsäure-Na-salz krist. reinst L-Lactic acid-Na-salt cryst. research grade	Kat.Nr./Cat.No. :	29763
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_3H_5O_3 \cdot Na$
Molekulargewicht molecular weight		112.1
Aussehen appearance		weißes körniges hygroskopisches Pulver white grained hygroscopic powder
Gehalt (%) assay	Titration	min. 99.0
Löslichkeit solubility	1 % in H ₂ O	klar und farblos clear and colorless
[α]²⁰_D	1 % in H ₂ O	-12.0° ± 1.0°
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.