

---

<b>1-Naphthylacetat p. A.</b> 1-Naphthyl acetate analytical grade	<b>Kat.Nr./Cat.No. :</b>	<b>30040</b>
--	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>12</sub> H <sub>10</sub> O <sub>2</sub>
<b>Molekulargewicht</b> molecular weight		186.21
<b>Aussehen</b> appearance		weiße bis fast weiße Kristalle white to off-white crystals
<b>Gehalt (%)</b> assay	HPLC	min. 99.0
<b>Schmelzpunkt (°C)</b> melting point		42 – 46
<b>A 1cm/ 280 nm</b>	0.001% in Methanol	min. 0.3
<b>Identität</b> identity	IR	entspricht corresponds
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.