
Nycodenz® 60% (w/v) Lösung in Wasser Nycodenz® 60% (w/v) solution in water	Kat.Nr./Cat.No. :	31000
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{19}H_{26}I_3N_3O_9$
Molekulargewicht molecular weight		821.1
Aussehen appearance		klare, farblose Lösung clear, colorless solution
Dichte (g/ml) density	20 °C	1.310 ± 0.002
Brechungsindex refractive index	20 °C	1.4273 ± 0.0003
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.