

Proteasome Inhibitor MG-132

Cat.No. :

33766

Parameter	Method	Specification
Amino acid-sequence		Z-Leu-Leu-Leu-CHO
Formula		C ₂₆ H ₄₁ N ₃ O ₅
Molecular weight		475.62
Appearance		white lyophilisate
Assay (%)	HPLC	min. 90
Solubility		soluble in DMSO
Advise		do not use an aqueous buffer substitute for DMSO because this compound may racemize soon
Storage / Stability	powder solution	at -15 to -25 °C in the dark for more than one year at -15 to -25 °C for more than one month if it is protected from moisture and kept in the dark

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.