

PRODUCT INFORMATION

GermDecon

Cat. No.: 34207

PRODUCT DESCRIPTION

Disinfectant for surfaces and laboratory instruments, efficient against a wide spectrum of microorganisms.

Application: GermDecon is for the use on any laboratory surfaces and equipment to avoid cross-contaminations or infections, especially in labs working with biological samples.

Decontamination of surfaces and laboratory instruments:

- Spray GermDecon directly on the surface/instrument to be disinfected until it is moistened completely.
 - Let the surface/instrument dry completely.
 - **IMPORTANT:** GermDecon does not leave any residues after drying.
-

Features

- Ready-to-use
 - Isopropyl alcohol based
 - Non-corrosive
 - Non-carcinogenic
-

Storage Store at + 15 °C to + 30 °C.

Ver. 07/15