

Serva Violett 17 reinst
Serva violet 17 research grade

Kat.Nr./Cat.No.: 35072

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₄₁ H ₄₄ N ₃ O ₆ S ₂ ·Na
Molekulargewicht molecular weight		761.9
Aussehen appearance		blauschwarzes Pulver blue-black powder
λ_{max} (nm)	0.001 % in H ₂ O	544 – 550
A 1cm/λ_{max}	0.001 % in H ₂ O	min. 0.8
Wasser (%) water	KF	max. 10.0
TLC	Kieselgel/silicagel 60 F ₂₅₄ n-Butanol/ n-butanol Ethanol/ ethanol Eisessig/ glacial acetic acid Wasser/ water	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.