

---

<b>3,3',5,5'-Tetramethylbenzidin</b> reinst 3,3',5,5'-Tetramethylbenzidine (TMB) research grade	<b>Kat.Nr./Cat.No. :</b>	<b>35926</b>
--	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>16</sub> H <sub>20</sub> N <sub>2</sub>
<b>Molekulargewicht</b> molecular weight		240.4
<b>Aussehen</b> appearance		hell elfenbeinfarbenes Pulver light ivory-colored powder
<b>Gehalt (%)</b> assay	Titration	98.0 – 102.0
<b>Schmelzpunkt (°C)</b> melting point		168 – 170
<b>Trocknungsverlust (%)</b> loss on drying		max. 0.5
<b>Identität</b> identity	IR	entspricht corresponds
<b>Eisen (ppm)</b> iron	ICP-OES	
<b>Lagerung (°C)</b> storage		+15 bis / to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.